



Narrow/Pediatric Socket

Ideal for children, teens and adults with narrow limbs

- New sleek lightweight socket
- Fit immediately without casts, scans or molds
- Updated closure system using high strength cord
- Comfort & stability with customizable, padded liner
- High strength, light weight socket, similar to aluminum substitutes
- Fit and aligned in a single visit with hand tools, ideal for any location
- Used as a preparatory, definitive or secondary device
- Waterproof and durable for any outdoor activity.
- Compatible with many commercially available feet, liners and components

Specifications

	Mobility Grade	K0 - K3
	Side	Left (L), Right (R)
	Distal Circumference (cm)	18 - 24 cm skin 20 - 26 cm with liner
	Base of Patella to End Length (cm)	Medium: Up to 20 cm Tall: Over 20 - 24 cm
	Product weight (Includes lock and pyramid)	1.4 lbs/ 635 grams
	Weight limit	Less than 110 kg / 220 lbs
	Closure System	Spool
	Water Capability	Waterproof
	Warranty	1 year
	Device Models	TTAC-N-M, TTAC-N-T

Prosthesis Kit



- Socket with cup/shuttle lock/ pyramid
- Neoprene liner (inside socket)
- Locking Pins 1.5"
- Loctite 242
- Two pull cords w/ cord (attached to socket)
- Padding kit (2 crescents, 2 rectangles, 1 large posterior pad)
- Cosmetic cover
- Sock
- Extras: pull tab, 4 feet of cord, 2 cord loops, 1 cleat and screws
- Tube clamp
- Pylon

Call Joshua Mullins at 301-821-5414 or email jmullins@ifitprosthetics.com for more information.

iFITProsthetics.com | Follow us:

iFIT Prosthetics is a Veteran-Owned Small Business (VOSB)

Patent #8470050, #8491667, #8845755, #10398577, #10806608, #11382775, #11766343. ©2023 iFIT Prosthetics, LLC. All rights reserved. Products mentioned herein are trademarks or registered trademarks of iFIT Prosthetics, LLC.



The Right Fit, Right Now™