














- Immediate fit with no casts, scans or molds—customizable to fit patient needs
- Comfort & stability with customizable, padded liner
- High strength, light weight socket, similar to aluminum substitutes
- Safety locking buckle system allows the user to adjust the socket size anytime
- Fit and aligned in a single visit with hand tools, ideal for any location
- Used as a preparatory, definitive or secondary device
- Perfect solution for patients with kidney or heart disease that experience daily volume changes
- Compatible with many commercially available feet, liners and components

Specifications

	Mobility Grade	K0 - K3
	Side	Left (L), Right (R)
	Distal Circumference (cm)	24 - 34
	Base of Patella to End Length (cm)	20 - 26
	Product weight	1520 g / 3.35 lbs
	Weight limit	118 kg / 260 lbs
	Buckles (Can be worn on either side)	Right or Left
	Closure System	Cable or Spool
	Water Capability	Waterproof
	Warranty	1 year
	Device Models	TT200Ri / TT200Li

Prosthesis Kit



- Socket
- Shuttle lock
- Pyramid adapter
- 2 Pins (2" and 2.5")
- Neoprene inner liner
- Padding kit
- Cables
- Power strap
- Pylon with bonded adapter
- Tube clamp
- Spacers

Call Joshua Mullins at 301-821-5414 or email jmullins@ifitprosthetics.com for more information.

iFITProsthetics.com | Follow us:   

 iFIT Prosthetics is a Veteran-Owned Small Business (VOSB)

Patent #8470050, #8491667, #8845755, #10398577, #10806608, #11382775. ©2023 iFIT Prosthetics, LLC. All rights reserved. Products mentioned herein are trademarks or registered trademarks of iFIT Prosthetics, LLC.



The Right Fit, Right Now™